



Responsible, Efficient, Reliable
PanSW Med Corp
www.panswmed.com

Pansw EKG cable

- 1, Various products compatible with many marks and models;
- 2, effective customers support for compatibilities and after sales service;
- 3, Various connector styles available: banana, clip, DIN, snap, etc;

DIN 3.0 connector



banana 4.0 end



Snap end



Clips end



R

N

F

L

C1

C2

C3

C4

C5

C6

IEC color coding

RA

RL

LL

LA

V1

V2

V3

V4

V5


V6

AHA color coding

Note: (within this catalogue)
1, integrated means complete cable from machine to electrode: trunk cable means cable from machine to the yoke; leadwire means the wires part from yoke to the electrodes;
2, EKG cable: means 10 leads patient cable for ECG machine or defibrillator;

Schiller

Part No	OEM No	Description		Picture	Compatibility
---------	--------	-------------	--	---------	---------------

EKG-10001A		10 leads; integrated; IEC, DB15M- 	<i>Schiller ECG: AT-102, CS-200, Cardiovit AT-10 plus, AT-104 PC, AT-110</i> <i>Schiller: AT1 ,AT2, AT3, AT5, AT6, AT10, AT60,AT104 PC,CS100,CS6, CS 6/12,AT4, AT104, AT101, AT102; Avionics(Del Mar): 910 / 920 / 930/ 940 / 942 / 950 / 960 / 970 / 980</i> <i>Custo-med: Custo-norm</i> <i>Cardioline: E1, E1R, E;Delta 3; Delta 300; C3</i> <i>Cardiorapid: K11, K111, K112, K131, K31, K33, K36</i> <i>Cardiotek: EP Tracer</i> <i>Dr Lee: ECG-310A; ECG-310B;</i> <i>Del mar: 910,920,930,940,942,950, 960,970,980</i> <i>Elettronica Trentina/Carditte: Autoruler 12/1; 12/3; Microruler 12/1; Microruler 12/3; Microtel: Excel 101; Excel 103; Excel 106;</i> <i>Ergoline: Autoruler, Autoscript 6/12</i> <i>Cardiette, EK 3003 / 3012, Ergoline</i> <i>Esaote Biomedica: EKG P120</i> <i>Fidelity: LP 3000</i> <i>Innomed: Heart Mirror1, Heart Mirror3; Heart Mirror3D; Heart Screen 80G; Heart Screen 112C; Heart Screen 112D;</i> <i>Kontron: EKG P80</i> <i>Progetti: EPG Smarty; EPG Pocket; EPG Junior; EPG View; EPG 6 View; EPG Plus; EPG 612; EPG 13</i> <i>Reynolds: Cardiotrack 12</i> <i>Schiller:AT3. AT5. AT6.AT10. AT60.</i> 
EKG-10001B		10 leads; integrated; AHA, DB15M- Di 	
EKG-10001C	AT-1	10 leads; integrated; IEC,DB15M- Banana 	
EKG-10001D		10 leads; integrated; AHA,DB15M- 	
EKG-10001E		10 leads; integrated; IEC,DB15M- 	
EKG-10001F		10 leads; integrated; AHA,DB15M- Clip end 	
EKG-10001G		10 leads; integrated; IEC,DB15M- snap e 	
EKG-10001H		10 leads; integrated; AHA,DB15M- 	
EKG-10001I		10 leads; ecg trunk cable, 10 leads, IEC; use with EKG-10045. Etc	
EKG-10001J		10 leads; ecg trunk cable, 10 leads, AHA color; use with EKG-10045, etc.	

HP Philips EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10040A		10 leads; integrated; IEC, DB15M- 		

EKG-10040B		10 leads; integrated; AHA, DB15M-			M1770A(Pagewriter 200i-300i); M171A (Pagewriter 200); M1772A (Pagewriter 100);	
EKG-10040C		10 leads; integrated; IEC,DB15M-				
EKG-10040D		10 leads; integrated; AHA,DB15M-				
EKG-10040E		10 leads; integrated; IEC,DB15M-				
EKG-10040F		10 leads; integrated; AHA,DB15M- Clip e				
EKG-10040G		10 leads; integrated; IEC,DB15M- snap e				
EKG-10040H		10 leads; integrated; AHA,DB15M- sna				
EKG-10040I		10 leads; ecg trunk cable that use along with GE ECG leads, IEC color				
EKG-10040J		10 leads; ecg trunk cable that use along with GE ECG leads, AHA color				

Philips EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10041A		10 leads lead set; IEC, Din3.0 		Compatible with: PageWriter Touch, 860284, PageWriter Trim Cardiographs,

EKG-10041B		10 leads lead set; AHA, Din3.0 			860286, 860288, 860290, 860297
EKG-10041C	98980 31291 91 (IEC)	10 leads lead set; IEC, Banana4.0 			
EKG-10041D	989803 129161 (AMMI/ AHA)	10 leads lead set; AHA, Banana4.0 			
EKG-10041E		10 leads lead set; IEC, clips 			
EKG-10041F		10 leads lead set; AHA, clips 			
EKG-10041G		10 leads lead set; AHA, snap 			
EKG-10041H		10 leads lead set; AHA, snap 			
EKG-10041I		Philips USB data cable			
EKG-10042A	98980 31291 71 (IEC color)	4 leads; limb leads; HP Pagewriter TRIM limb 4 lead set, banana 4.0 end; IEC			IEC 989803129171
EKG-10042B	98980 31291 41 (AHA color)	4 leads; limb leads; HP Pagewriter TRIM limb 4 lead set, banana 4.0 end; AHA/AAMI			AAMI 989803129141
EKG-10042C	98980 31291 81 (IEC color)	6 leads, chest leads; HP Pagewriter trim Chest 6 lead set, banana 4.0 end; IEC			IEC 989803129181

EKG-10042D	98980 31291 51 (AHA color)	6 leads, chest leads; HP Pagewriter trim Chest 6 lead set, banana 4.0 end; AHA/AAMI			AAMI 989803129151
------------	--	---	--	--	-------------------

Kenz

Part No	OEM No	Description	Connector	Compatibility
EKG-10003A	PC-104	10 leads; integrated; IEC, DB15M- 		
EKG-10003B		10 leads; integrated; AHA, DB15M- 		
EKG-10003C		10 leads; integrated; IEC, DB15M- 		
EKG-10003D		10 leads; integrated; AHA, DB15M- 		
EKG-10003E		10 leads; integrated; IEC, DB15M- 		
EKG-10003F		10 leads; integrated; AHA, DB15M- Clip end 		
EKG-10003G		10 leads; integrated; IEC, DB15M- snap end 		
EKG-10003H		10 leads; integrated; AHA, DB15M- snap end 		

*Kenz Delta
plus, delta
3plus, delta 60
plus, Cardioline
Delta 1- Delta 1
Plus
Delta 3-Delta 3
Plus Delta 30 Delta
60-Delta 60 Plus*

Cardioline; Kenz

Part No	OEM No	Description		Connector	Compatibility
EKG-10004A		10 leads; integrated; IEC, DB15M- 		 PanSW Med Corp www.panswmed.com	Cardioline 108/110,1203,1205; Suzuken kenz
EKG-10004B		10 leads; integrated; AHA, DB15M- 			
EKG-10004C	PC-109	10 leads; integrated; IEC,DB15M-Banana4.0 			
EKG-10004D		10 leads; integrated; AHA,DB15M-Banana4.0 			
EKG-10004E		10 leads; integrated; IEC,DB15M- 			
EKG-10004F		10 leads; integrated; AHA,DB15M-Clip end 			
EKG-10004G		10 leads; integrated; IEC,DB15M-snap end 			
EKG-10004H		10 leads; integrated; AHA,DB15M-snap end 			


Nihon Kohden

Part No	OEM No	Description		Connector	Compatibility
EKG-10005A	K079 (BJ-901D)	10 leads; integrated; IEC, DB15M-Din3.0 			

EKG-10005B	K019 (BA-901D)	10 leads; integrated; AHA, DB15M-Din3.0 		<i>Nihon Kohden ECG-9620, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9130, ECG-9132, DongJiang ECG-11A, ECG-11B, ECG-11D, ECG-32A, (II), ECG-120 For ECG-1150/1250/1350/2350</i>
EKG-10005C	K081 (BJ-902D)	10 leads; integrated; IEC, DB15M-Banana4.0 		
EKG-10005D		10 leads, integrated; AHA, DB15M-Banana4.0 		
EKG-10005E		10 leads, integrated; IEC, DB15M-Clip end 		
EKG-10005F		10 leads, integrated; AHA, DB15M-Clip end 		
EKG-10005G	K082 BJ-903D	10 leads; integrated; IEC, DB15M-snap 		
EKG-10005H		10 leads; integrated; AHA, DB15M-snap 		
EKG-10005I	BJ-904D; K098	5 leads integrated; IEC, din3.0		for ECG-1950
EKG-10005J		5 leads integrated; AHA, din3.0		
EKG-10005K		5 leads integrated; IEC, banana 		

EKG-10005L				
5 leads integrated; IEC, banana				

Nihon Kohden

Part No	OEM No	Description	Connector	Compatibility
EKG-10006A		10 leads; integrated; IEC, DB15M-Din3.0 		for ECG-1550
EKG-10006B		10 leads; integrated; AHA, DB15M-Din3.0 		
EKG-10006C	BJ-912D	10 leads; integrated; IEC, DB15M-Banana4.0 		
EKG-10006D		10 leads; integrated; AHA, DB15M-Banana4.0 		
EKG-10006E	BJ-913D	10 leads; integrated; IEC, DB15M-Clip end 		
EKG-10006F		10 leads; integrated; AHA, DB15M-Clip end 		
EKG-10006G		10 leads; integrated; IEC, DB15M-snap end 		
EKG-10006H		10 leads; integrated; AHA, DB15M-snap end 		

Shanghai Kohden

Part No	OEM No	Description	Connector	Compatibility
---------	--------	-------------	-----------	---------------

EKG-10007A		10 leads; integrated; IEC, DB15M- 		Shanghai Kohden:ECG- 6511,ECG- 6151,8110P/K,635 3
EKG-10007B		10 leads; integrated; AHA, DB15M- 		
EKG-10007C		10 leads; integrated; IEC,DB15M- 		
EKG-10007D		10 leads; integrated; AHA,DB15M- 		
EKG-10007E		10 leads; integrated; IEC,DB15M- 		
EKG-10007F		10 leads; integrated; AHA,DB15M- Clip e 		
EKG-10007G		10 leads; integrated; IEC,DB15M- s 		
EKG-10007H		10 leads; integrated; AHA,DB15M- snap 		

Burdick

Part No	OEM No	Description		Connector	Compatibility
EKG-10008A	EK-10	10 leads; integrated; IEC, DB15M-Din 			
EKG-10008B		10 leads; integrated; AHA, DB15M-Din3.0 			

EKG-10008C		10 leads; integrated; IEC,DB15M-Banana4.0		    	Burdick EK-10
EKG-10008D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10008E		10 leads; integrated; IEC,DB15M-			
EKG-10008F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10008G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10008H		10 leads; integrated; AHA,DB15M-snap end			



Burdick Quinton Eclipse EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10051A		10 leads; integrated; IEC, DB15M-Din3.0	   	 Burdick Quinton Eclipse 4 3 100 800
EKG-10051B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10051C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10051D		10 leads; integrated; AHA,DB15M-Banana4.0		

EKG-10051E		10 leads; integrated; IEC,DB15M-Clip end 			4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS; burdick 8500
EKG-10051F		10 leads; integrated; AHA,DB15M-Clip end 			
EKG-10051G		10 leads; integrated; IEC,DB15M-snap end 			
EKG-10051H		10 leads; integrated; AHA,DB15M-snap end 			

Fukuda

Part No	OEM No	Description	Connector	Compatibility
EKG-10009A	FX-101	10 leads; integrated; IEC, DB15M-Din3.0 		Fukuda FCP-7411, FX-7402, FX-7202, FX-2111,FX-2155, FX-3010, FX-7102
EKG-10009B		10 leads; integrated; AHA, DB15M-Din3.0 		
EKG-10009C		10 leads; integrated; IEC,DB15M-Banana4.0 		
EKG-10009D		10 leads; integrated; AHA,DB15M-Banana4.0 		
EKG-10009E		10 leads; integrated; IEC,DB15M-Clip end 		
EKG-10009F		10 leads; integrated; AHA,DB15M-Clip end 		

EKG-10009G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10009H		10 leads; integrated; AHA,DB15M-snap end			
Part No	OEM No	Description		Connector	Compatibility
EKG-10052A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10052B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10052C		10 leads; integrated; IEC,DB15M-Banana4.0			
EKG-10052D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10052E		10 leads; integrated; IEC,DB15M-Clip end			KP-309/KP-405/KP-500/ C100 C110 C120 C300 C310 C550 C600 C600W C610 C610W C625 a2000 series; a1000
EKG-10052F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10052G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10052H		10 leads; integrated; AHA,DB15M-snap end			
Part No	OEM No	Description		Connector	Compatibility


EKG-10053A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10053B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10053C		10 leads; integrated; IEC, DB15M-Banana4.0			
EKG-10053D		10 leads; integrated; AHA, DB15M-Banana4.0			
EKG-10053E		10 leads; integrated; IEC, DB15M-Clip end			
EKG-10053F		10 leads; integrated; AHA, DB15M-Clip end			
EKG-10053G		10 leads; integrated; IEC, DB15M-snap end			
EKG-10053H		10 leads; integrated; AHA, DB15M-snap end			




KP-309D/KP-405D/KP-500D/
C100 C110 C120
C300 C310 C550
C600 C600W C610
C610W C625
a2000 series;
a1000; (with 20k
built-in resistor for
defibrillation
protection)

GE-Marquette

Part No	OEM No	Description	Connector	Compatibility
EKG-10043A		10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10043B		10 leads; integrated; AHA, DB15M-Din3.0		

EKG-10043C	20298 90-001	10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10043D	20298 93-001	10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10043E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10043F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10043G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10043H		10 leads; integrated; AHA,DB15M-snap end		
EKG-10044A	22341 808	Trunk cable, IEC		
EKG-10044B		Trunk cable, AHA		
EKG-10045A		10 leads lead set; IEC, Din3.0		
EKG-10045B		10 leads lead set; AHA, Din3.0		
EKG-10045C		10 leads lead set; IEC, Banana4.0		



*Hellige:
MicroSmart, MAC
500/1100; MAC
1200 / 1200 ST*

EKG-10045D	10 leads lead set; AHA, Banana4.0		
EKG-10045E	10 leads lead set; IEC, clips		
EKG-10045F	10 leads lead set; AHA, clips		
EKG-10045G	10 leads lead set; AHA, snap		
EKG-10045H	10 leads; lead set; AHA, snap		
EKG-10046A	10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10046B	10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10046C	10 leads; integrated; IEC, DB15M-Banana4.0		
EKG-10046D	10 leads; integrated; AHA, DB15M-Banana4.0		
EKG-10046E	10 leads; integrated; IEC, DB15M-Clip end		
EKG-10046F	10 leads; integrated; AHA, DB15M-Clip end		

EKG-10046G		10 leads; integrated; IEC,DB15M-snap end			For GE Eagle, Solar, Dash , Tram MAC-Lab Cath Lab System, Datex-Ohmeda S/5 FM, Original P/N: (AAMI: E9002ZH, 416035-001) (IEC: E9002ZJ, 416035-002)
EKG-10046H		10 leads; integrated; AHA,DB15M-snap end			
EKG-10047A		10 leads; Trunk cable, IEC			
EKG-10047B		10 leads; Trunk cable, AHA			
EKG-10048A		10 leads lead set; IEC, Din3.0			
EKG-10048B		10 leads lead set; AHA, Din3.0			
EKG-10048C		10 leads lead set; IEC, Banana4.0			
EKG-10048D		10 leads lead set; AHA, Banana4.0			
EKG-10048E		10 leads lead set; IEC, clips		414556-004 Set/5, Grabber connector, R,L,F,N,C, 130 cm, IEC; 416467-004 74 cm grabber connector, C2-- C6, IEC 416467-002 130cm, C2-C6 IEC	
EKG-10048F		10 leads lead set; AHA, clips			
EKG-10048G		10 leads lead set; AHA, snap			

EKG-10048H		10 leads; lead set; AHA, snap		
------------	--	-------------------------------	---	--

NEC

Part No	OEM No	Description	Connector	Compatibility
EKG-10011A		10 leads; integrated; IEC, DB15M-Din3.0 		NEC:K-3321,3330,3370
EKG-10011B		10 leads; integrated; AHA, DB15M-Din3.0 		
EKG-10011C		10 leads; integrated; IEC, DB15M-Banana4.0 		
EKG-10011D		10 leads; integrated; AHA, DB15M-Banana4.0 		
EKG-10011E		10 leads; integrated; IEC, DB15M-Clip end 		
EKG-10011F		10 leads; integrated; AHA, DB15M-Clip end 		
EKG-10011G		10 leads; integrated; IEC, DB15M-snap end 		
EKG-10011H		10 leads; integrated; AHA, DB15M-snap end 		

Suction Electrode

Part No	OEM No	Description		Connector	Patient group
EKG-10012A	21714 401	Suction Electrode use with DIN 3.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult
EKG-10012B		Suction Electrode use with banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012C		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012D	20298 92- 001; Philips 40418 A	Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 26mm; metal electrode diameter: 21 mm;			adult size
EKG-10012E		Suction Electrode use with DIN 3.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012F		Suction Electrode use with banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012G		suction ECG electrode use with both DIN 3.0 and banana 4.0 connector; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size
EKG-10012H		Multi-Use Suction Electrode use for snap/banana/DIN/clip connectors; ball diameter: 23mm; metal electrode diameter: 17mm;			pediatric size

Limb ECG Electrode

Part No	OEM No	Description		Connector	Compatibility
---------	--------	-------------	--	-----------	---------------

EKG-10013A		Limb Electrode for DIN 3.0 connector; clamps size: 30x130mm;			for adult use
EKG-10013B		Limb Electrode for banana 4.0 connector; clamps size: 30x130mm;			for adult use
EKG-10013C		limb Electrode for DIN 3.0 connector and banana 4.0 connector; clamps size: 30x130mm;			for adult use
EKG-10013D	20298 91-001; Philips 40494 E	Multi-Use Limb Electrode; for snap/DIN/banana/clip connectors; universal use; clamps size: 30x130mm;		 	for adult use
EKG-10013E		Multi-Use Limb Electrode; for snap/DIN/banana/clip connectors; universal use; clamps size: 20x85mm			for pediatric use

SIEMENS

EKG-10015A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10015B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10015C		10 leads; integrated; IEC, DB15M-Banana4.0			
EKG-10015D		10 leads; integrated; AHA, DB15M-Banana4.0			
EKG-10015E		10 leads; integrated; IEC, DB15M-Clip end			Cardiostat 31

EKG-10015F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10015G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10015H		10 leads; integrated; AHA,DB15M-snap end		

Bionet

EKG-10016A		10 leads; integrated; IEC, DB15M-Din3.0			<i>Cardiocare 2000, Cardiotouch, Medical Econet Cardio M y Cardio M Plus</i>
EKG-10016B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10016C		10 leads; integrated; IEC,DB15M-Banana4.0			
EKG-10016D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10016E		10 leads; integrated; IEC,DB15M-Clip end			
EKG-10016F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10016G		10 leads; integrated; IEC,DB15M-snap end			

EKG-10016H		10 leads; integrated; AHA,DB15M-snap end		
------------	--	--	---	--

MORTARA

EKG-10018A		10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10018B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10018C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10018D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10018E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10018F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10018G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10018H	9293-022-50	10 leads; integrated; AHA,DB15M-snap end		

MORTARA
INSTRUMENT:ELI
100 MORTARA
INSTRUMENT:ELI
100 MORTARA
INSTRUMENT:ELI
200 MORTARA
INSTRUMENT:ELI
50

Mortara Eli 250 cables


Part No	OEM No	Description	Connector	Compatibility
EKG-10054A		10 leads; integrated; IEC, DB15M-Din3.0		

EKG-10054B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10054C		10 leads; integrated; IEC,DB15M-Banana4.0			
EKG-10054D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10054E		10 leads; integrated; IEC,DB15M-Clip end			
EKG-10054F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10054G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10054H		10 leads; integrated; AHA,DB15M-snap end			

MORTARA
INSTRUMENT: ELI
10/ELI150/ELI250/
ELI210/ELI350

Cardioline; Cardiette

Part No	OEM No	Description		Connector	Compatibility
EKG-10019A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10019B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10019C		10 leads; integrated; IEC,DB15M-Banana4.0			

EKG-10019D		10 leads; integrated; AHA,DB15M-Banana4.0		AR-600, AR-1200, part number 63050025
EKG-10019E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10019F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10019G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10019H		10 leads; integrated; AHA,DB15M-snap end		

Edan

Part No	OEM No	Description	Connector	Compatibility
EKG-10020A		10 leads; integrated; IEC, DB15M-Din3.0		SE-12 Express SE-3 SE-601A
EKG-10020B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10020C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10020D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10020E		10 leads; integrated; IEC,DB15M-Clip end		

EKG-10020F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10020G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10020H		10 leads; integrated; AHA,DB15M-snap end		

Hellige

Part No	OEM No	Description	Picture	Comaptible
EKG-10021A		10 leads; integrated; IEC, DB15M-Din3.0		 Bosch:603D,603ER,603S,606 Cambridge: CM 6000, CM 6400 Custo-med: Custo-card DEGO: Cardiomed 3 Hellige: EK36,EK41,EK43,EK53, EK 56, EK 403, EK 413, EK 512, EK 512E,EK 512 P, Cardiosmart, CardioSys, Corna Siemens: Cardiostat1/3T/703/710,Cardina C, Cardirex 3T/4/6T/62,Cathcor, Megacard ElemaA/B,R,R/E, Mingodoc,Mingocard 3/4/7,Mingograph 4/62/81 Mingograph 82/740,Mingolog7,Mingophon3/7/310/410/420/710/720,Miracard 1, Sirecard440/440S/460/460S/740Siereg, Sirecard 803,Sirecust 701; SICARD 440
EKG-10021B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10021C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10021D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10021E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10021F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10021G		10 leads; integrated; IEC,DB15M-snap end		

EKG-10021H		10 leads; integrated; AHA,DB15M-snap end		
------------	--	--	---	--

Medtronic/Physio Control

Part No	OEM No	Description	Picture	Compatible
EKG-10022A		10 leads; integrated; IEC, DB15M-Din3.0		<i>Lifepak 12</i>
EKG-10022B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10022C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10022D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10022E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10022F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10022G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10022H		10 leads; integrated; AHA,DB15M-snap end		

Dixtal EKG cable

Part No	OEM No	Description	Picture	Compatible
---------	--------	-------------	---------	------------

EKG-10023A		10 leads; integrated; IEC, DB15M-Din3.0				EP-3; DX-203
EKG-10023B		10 leads; integrated; AHA, DB15M-Din3.0				
EKG-10023C		10 leads; integrated; IEC, DB15M-Banana4.0				
EKG-10023D		10 leads; integrated; AHA, DB15M-Banana4.0				
EKG-10023E		10 leads; integrated; IEC, DB15M-Clip end				
EKG-10023F		10 leads; integrated; AHA, DB15M-Clip end				
EKG-10023G		10 leads; integrated; IEC, DB15M-snap end				
EKG-10023H		10 leads; integrated; AHA, DB15M-snap end				

Dixtal EP-12 cables					
Part No	OEM No	Description		Connector	Compatibility
EKG-10055A		10 leads; integrated; IEC, DB15M-Din3.0			EP-12
EKG-10055B		10 leads; integrated; AHA, DB15M-Din3.0			

EKG-10055C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10055D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10055E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10055F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10055G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10055H		10 leads; integrated; AHA,DB15M-snap end		

Welch Allyn EKG cable

Part No	OEM No	Description	Picture	Comaptible
EKG-10024A		10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10024B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10024C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10024D		10 leads; integrated; AHA,DB15M-Banana4.0		

EKG-10024E	10 leads; integrated; IEC,DB15M-Clip end		  	Welch Allyn Cardioperfect MD
EKG-10024F	10 leads; integrated; AHA,DB15M-Clip end			
EKG-10024G	10 leads; integrated; IEC,DB15M-snap end			
EKG-10024H	10 leads; integrated; AHA,DB15M-snap end			
EKG-10056A	10 leads; integrated; IEC, Lemo 14P-Din3.0		    	Welch Allyn Cardioperfect Pro
EKG-10056B	10 leads; integrated; AHA, Lemo 14P-Din3.0			
EKG-10056C	10 leads; integrated; IEC,Lemo 14P-Banana4.0			
EKG-10056D	10 leads; integrated; AHA,Lemo 14P-Banana4.0			
EKG-10056E	10 leads; integrated; IEC,Lemo 14P-Clip end			
EKG-10056F	10 leads; integrated; AHA,Lemo 14P-Clip end			
EKG-10056G	10 leads; integrated; IEC,Lemo 14P-snap end			




EKG-10056H		10 leads; integrated; AHA,Lemo 14P-snap end			
EKG-10057A		10 leads lead set; IEC, Din3.0			
EKG-10057B		10 leads lead set; AHA, Din3.0			
EKG-10057C		10 leads lead set; IEC, Banana4.0			
EKG-10057D	401129	10 leads lead set; AHA, Banana4.0			
EKG-10057E		10 leads lead set; IEC, clips			
EKG-10057F		10 leads lead set; AHA, clips			
EKG-10057G		10 leads lead set; AHA, snap			
EKG-10057H		10 leads lead set; AHA, snap			


welch allyn CP100
CP200 ECG cable
patient leads

Nihon Kohden

Part No	OEM No	Description	Connector	Compatibility
EKG-10025A		10 leads; integrated; IEC, DB15M-Din3.0		

EKG-10025B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10025C		10 leads; integrated; IEC,DB15M-Banana4.0			
EKG-10025D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10025E	BJ900P	10 leads; integrated; IEC,DB15M-Clip end			
EKG-10025F	BJ900 PA	10 leads; integrated; AHA,DB15M-Clip end			
EKG-10025G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10025H		10 leads; integrated; AHA,DB15M-snap end			

Philips ECG cable, new ten leads

Part No	OEM No	Description	Connector	Compatibility
EKG-10026A		10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10026B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10026C		10 leads; integrated; IEC,DB15M-Banana4.0		

EKG-10026D		10 leads; integrated; AHA,DB15M-Banana4.			
EKG-10026E		10 leads; integrated; IEC,DB15M-Clip end			
EKG-10026F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10026G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10026H		10 leads; integrated; AHA,DB15M-snap end			
EKG-10026I	9898 0314 7691 ; M1 949A	10 leads; ECG trunk cable M1949A; M1663A			HeartStart MRx: M3535A, M3536A

Ecafix Funbec 12S

Part No	OEM No	Description	Connector	Compatibility
EKG-10027A		10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10027B		10 leads; integrated; AHA, DB15M-Din3.0		

EKG-10027C		10 leads; integrated; IEC,DB15M-Banana4.0		Funbec 12S
EKG-10027D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10027E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10027F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10027G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10025H		10 leads; integrated; AHA,DB15M-snap end		

MICROMED

Part No	OEM No	Description	Connector	Compatibility
EKG-10028A		11 leads; integrated; IEC, DB15M-Din3.0		
EKG-10028B		11 leads; integrated; AHA, DB15M-Din3.0		
EKG-10028C		11 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10028D		11 leads; integrated; AHA,DB15M-Banana4.0		

EKG-10028E		11 leads; integrated; IEC,DB15M-Clip end		
EKG-10028F		11 leads; integrated; AHA,DB15M-Clip end		
EKG-10028G		11 leads; integrated; IEC,DB15M-snap end		
EKG-10028H		11 leads; integrated; AHA,DB15M-snap end		

Carewell

Part No	OEM No	Description	Connector	Compatibility
EKG-10029A		10 leads; integrated; IEC, DB15M-Din3.0		Carewell All, 1103 1106...
EKG-10029B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10029C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10029D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10029E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10029F		10 leads; integrated; AHA,DB15M-Clip end		

EKG-10029G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10029H		10 leads; integrated; AHA,DB15M-snap end		

BURDICK quinton EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10030A		10 leads; integrated; IEC, DB15M-Din3.0		<i>BURDICK QUINTON Q STRESS</i>
EKG-10030B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10030C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10030D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10030E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10030F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10030G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10030H		10 leads; integrated; AHA,DB15M-snap end		

Nelch Allyn PIC50 EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10031A		10 leads; integrated; IEC, DB15M-Din3.0 	 PanSW Med Corp www.panswmed.com	welch allyn PIC 50 ten leads cable
EKG-10031B		10 leads; integrated; AHA, DB15M-Din3.0 		
EKG-10031C		10 leads; integrated; IEC, DB15M-Banana4.0 		
EKG-10031D		10 leads; integrated; AHA, DB15M-Banana4.0 		
EKG-10031E		10 leads; integrated; IEC, DB15M-Clip end 		
EKG-10031F		10 leads; integrated; AHA, DB15M-Clip end 		
EKG-10031G		10 leads; integrated; IEC, DB15M-snap end 		
EKG-10031H		10 leads; integrated; AHA, DB15M-snap end 		

mortara EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10032A		10 leads; integrated; IEC, DB15M-Din3.0 	ROUND 12 pins	


EKG-10032B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10032C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10032D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10032E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10032F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10032G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10032H		10 leads; integrated; AHA,DB15M-snap end		

Mindray EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10033A		10 leads; integrated; IEC, DB15M-Din3.0		
EKG-10033B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10033C		10 leads; integrated; IEC,DB15M-Banana4.0		

EKG-10033D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10033E		10 leads; integrated; IEC,DB15M-Clip end			
EKG-10033F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10033G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10033H		10 leads; integrated; AHA,DB15M-snap end			
EKG-10033I		10 leads; AHA/IEC, TRUNK CABLE			

Philips TC30 TC50 ECG leads

Part No	OEM No	Description		Connector	Compatibility
EKG-10034A		10 leads lead set; IEC, Din3.0			
EKG-10034B		10 leads lead set; AHA, Din3.0			
EKG-10034C	98980 31516 51	10 leads lead set; IEC, Banana4.0			Philips EKG


EKG-10034D	989803 151631	10 leads lead set; AHA, Banana4.0 		philips EKG machine TC30, TC50, TC70
EKG-10034E		10 leads lead set; IEC, clips 		
EKG-10034F		10 leads lead set; AHA, clips 		
EKG-10034G		10 leads lead set; AHA, snap 		
EKG-10034H		10 leads lead set; AHA, snap 		

Mortara EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10035A		10 leads lead set; IEC, Din3.0 		ELI150C ELI230 ELI250C ELI350 wireless acquisition module ECG cable
EKG-10035B		10 leads lead set; AHA, Din3.0 		
EKG-10035C		10 leads lead set; IEC, Banana4.0 		
EKG-10035D		10 leads lead set; AHA, Banana4.0 		
EKG-10035E		10 leads lead set; IEC, clips 		

EKG-10035F		10 leads lead set; AHA, clips 			
EKG-10035G		10 leads lead set; AHA, snap 			
EKG-10035H		10 leads lead set; AHA, snap 			

Zoll EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10036A		10 leads; integrated; IEC, DB15M-Din3.0 		Zoll
EKG-10036B		10 leads; integrated; AHA, DB15M-Din3.0 		
EKG-10036C		10 leads; integrated; IEC, DB15M-Banana4.0 		
EKG-10036D		10 leads; integrated; AHA, DB15M-Banana4.0 		
EKG-10036E		10 leads; integrated; IEC, DB15M-Clip end 		
EKG-10036F		10 leads; integrated; AHA, DB15M-Clip end 		
EKG-10036G		10 leads; integrated; IEC, DB15M-snap end 		


EKG-10036H		10 leads; integrated; AHA,DB15M-snap end		
------------	--	--	---	--

Spacelabs EKG cable

Part No	OEM No	Description	Connector	Compatibility
EKG-10037A		10 leads; integrated; IEC, DB15M-Din3.0		Spacelabs
EKG-10037B		10 leads; integrated; AHA, DB15M-Din3.0		
EKG-10037C		10 leads; integrated; IEC,DB15M-Banana4.0		
EKG-10037D		10 leads; integrated; AHA,DB15M-Banana4.0		
EKG-10037E		10 leads; integrated; IEC,DB15M-Clip end		
EKG-10037F		10 leads; integrated; AHA,DB15M-Clip end		
EKG-10037G		10 leads; integrated; IEC,DB15M-snap end		
EKG-10037H		10 leads; integrated; AHA,DB15M-snap end		

Biocare EKG cable

Part No	OEM No	Description	Connector	Compatibility
---------	--------	-------------	-----------	---------------

EKG-10038A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10038B		10 leads; integrated; AHA, DB15M-Din3.0			
EKG-10038C		10 leads; integrated; IEC, DB15M-Banana4.0			
EKG-10038D		10 leads; integrated; AHA, DB15M-Banana4.0			
EKG-10038E		10 leads; integrated; IEC, DB15M-Clip end			
EKG-10038F		10 leads; integrated; AHA, DB15M-Clip end			
EKG-10038G		10 leads; integrated; IEC, DB15M-snap end			
EKG-10038H		10 leads; integrated; AHA, DB15M-snap end			

Bexen Cardio Reanibex 800 EKG cable

Part No	OEM No	Description		Connector	Compatibility
EKG-10039A		10 leads; integrated; IEC, DB15M-Din3.0			
EKG-10039B		10 leads; integrated; AHA, DB15M-Din3.0			

EKG-10039C		10 leads; integrated; IEC,DB15M-Banana4.0			
EKG-10039D		10 leads; integrated; AHA,DB15M-Banana4.0			
EKG-10039E		10 leads; integrated; IEC,DB15M-Clip end			
EKG-10039F		10 leads; integrated; AHA,DB15M-Clip end			
EKG-10039G		10 leads; integrated; IEC,DB15M-snap end			
EKG-10039H		10 leads; integrated; AHA,DB15M-snap end			

Philips EKG cable

Part No	OEM No	Description		Connector	Compatibility
EKG-10040A		integrated; IEC, DB15M-Din3.0			M1770A(Pagewriter 200i-300i); M171A (Pagewriter 200); M1772A (Pagewriter 100);
EKG-10040B		integrated; AHA, DB15M-Din3.0			
EKG-10040C		integrated; IEC, DB15M-Banana4.0			
EKG-10040D		integrated; AHA, DB15M-Banana4.0			
EKG-10040E		integrated; IEC, 10 leads clip end			
EKG-10040F		integrated; AHA, DB15M-clip			
EKG-10040G		integrated; IEC, DB15, snap			
EKG-10040H		integrated; AHA, DB15M-snap			
EKG-10040I		ecg trunk cable that use along with GE ECG leads, IEC color			

EKG-10040J		ecg trunk cable that use along with GE ECG leads, AHA color		
------------	--	---	--	--

Philips EKG cable

Part No	OEM No	Description		Connector	Compatibility
EKG-10041A		10 leads lead set; IEC, Din3.0		 	Compatible with: PageWriter Touch, 860284, PageWriter Trim Cardiographs, 860286, 860288, 860290, 860297
EKG-10041B		10 leads lead set; AHA, Din3.0			
EKG-10041C		10 leads lead set; IEC, Banana4.0			
EKG-10041D		10 leads lead set; AHA, Banana4.0			
EKG-10041E		10 leads lead set; IEC, clips			
EKG-10041F		10 leads lead set; AHA, clips			
EKG-10041G		10 leads lead set; AHA, snap			
EKG-10041H		10 leads lead set; AHA, snap			
EKG-10041I		Philips USB data cable		 	
EKG-10042A		HP Pagewriter TRIM limb 4 lead set, banana 4.0 end; IEC			IEC 989803129171
EKG-10042B		HP Pagewriter TRIM limb 4 lead set, banana 4.0 end; AHA/AAMI			AAMI 989803129141
EKG-10042C		HP Pagewriter trim Chest 6 lead set, banana 4.0 end; IEC			IEC 989803129181
EKG-10042D		HP Pagewriter trim Chest 6 lead set, banana 4.0 end; AHA/AAMI			AAMI 989803129151
Part No	OEM No	Description		Connector	Compatibility
EKG-10043A		integrated; IEC, DB15M-Din3.0		  	
EKG-10043B		integrated; AHA, DB15M-Din3.0			
EKG-10043C	20298 90-001	integrated; IEC, DB15M-Banana4.0			
EKG-10043D	20298 93-001	integrated; AHA, DB15M-Banana4.0			
EKG-10043E		integrated; IEC, 10 leads clip end			
EKG-10043F		integrated; AHA, DB15M-clip			
EKG-10043G		integrated; IEC, DB15, snap			
EKG-10043H		integrated; AHA, DB15M-snap			
EKG-10044A	22341	Trunk cable, IEC			Hellige: MicroSmart, MAC 500/1100; MAC 1200, 1200 CT

EKG-10044B	808	Trunk cable, AHA		1200 / 1200 SI
EKG-10045A		10 leads lead set; IEC, Din3.0	 	
EKG-10045B		10 leads lead set; AHA, Din3.0		
EKG-10045C		10 leads lead set; IEC, Banana4.0		
EKG-10045D		10 leads lead set; AHA, Banana4.0		
EKG-10045E		10 leads lead set; IEC, clips		
EKG-10045F		10 leads lead set; AHA, clips		
EKG-10045G		10 leads lead set; AHA, snap		
EKG-10045H		10 leads lead set; AHA, snap		
EKG-10046A		integrated; IEC, DB15M-Din3.0		
EKG-10046B		integrated; AHA, DB15M-Din3.0	 	
EKG-10046C		integrated; IEC, DB15M-Banana4.0		
EKG-10046D		integrated; AHA, DB15M-Banana4.0		
EKG-10046E		integrated; IEC, 10 leads clip end		
EKG-10046F		integrated; AHA, DB15M-clip		
EKG-10046G		integrated; IEC, DB15, snap		
EKG-10046H		integrated; AHA, DB15M-snap		
EKG-10047A		Trunk cable, IEC		For GE Eagle, Solar, Dash , Tram MAC-Lab Cath Lab System, Datex-Ohmeda S/5 FM, Original P/N: (AAMI: E9002ZH, 416035-001) (IEC: E9002ZJ, 416035-002)
EKG-10047B		Trunk cable, AHA		
EKG-10048A		10 leads lead set; IEC, Din3.0		
EKG-10048B		10 leads lead set; AHA, Din3.0		
EKG-10048C		10 leads lead set; IEC, Banana4.0		
EKG-10048D		10 leads lead set; AHA, Banana4.0		
EKG-10048E		10 leads lead set; IEC, clips	414556-004 Set/5, Grabber connector, R,L,F,N,C, 130 cm, IEC; 416467-004 74 cm grabber connector, C2-- C6, IEC 416467-002 130cm, C2-C6 IEC	
EKG-10048F		10 leads lead set; AHA, clips		
EKG-10048G		10 leads lead set; AHA, snap		
EKG-10048H		10 leads lead set; AHA, snap		

Medtronic EKG cable

Part No	OEM No	Description		Connector	Compatibility
EKG-10049		EKG trunk cable for medtronic; it works for both IEC & AHA color; cable length 2.2 meters			Lifepak 12
					

Medtronic EKG cable

Part No	OEM No	Description		Connector	Compatibility
EKG-10050A		6 leads set IEC, Din3.0			
EKG-10050B		6 leads set AHA, Din3.0			
EKG-10050C		6 leads set IEC, Banana4.0			
EKG-10050D		6 leads set AHA, Banana4.0			
EKG-10050E		6 leads set IEC, clip end			
EKG-10050F		6 leads set AHA, clip			
EKG-10050G		6 leads set IEC, snap			
EKG-10050H		6 leads set AHA, snap			

Burdick Quinton Eclipse EKG cable

Part No	OEM No	Description		Connector	Compatibility
EKG-10051A		10 leads integrated; IEC, DB15M-Din3.0			
EKG-10051B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10051C		10 leads integrated; IEC, DB15M-Banana4.0			

EKG-10051D		10 leads integrated; AHA,DB15M-Banana4.0			<i>Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS; burdick 8500</i>
EKG-10051E		10 leads integrated; IEC,DB15M-Clip end			
EKG-10051F		10 leads integrated; AHA,DB15M-Clip end			
EKG-10051G		10 leads integrated; IEC,DB15M-snap end			
EKG-10051H		10 leads integrated; AHA,DB15M-snap end			
Part No	OEM No	Description		Connector	Compatibility
EKG-10052A		10 leads integrated; IEC, DB15M-Din3.0			<i>KP-309/KP-405/KP-500/ C100 C110 C120 C300 C310 C550 C600 C600W C610 C610W C625 a2000 series; a1000</i>
EKG-10052B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10052C		10 leads integrated; IEC,DB15M-Banana4.0			
EKG-10052D		10 leads integrated; AHA,DB15M-Banana4.0			
EKG-10052E		10 leads integrated; IEC,DB15M-Clip end			
EKG-10052F		10 leads integrated; AHA,DB15M-Clip end			

EKG-10052G		10 leads integrated; IEC,DB15M-snap end 			
EKG-10052H		10 leads integrated; AHA,DB15M-snap end 			
Part No	OEM No	Description		Connector	Compatibility
EKG-10053A		10 leads integrated; IEC, DB15M-Din 			KP-309D/KP-405D/KP-500D/ C100 C110 C120 C300 C310 C550 C600 C600W C610 C610W C625 a2000 series; a1000; with 20k built-in resistor for defibrillation protection)
EKG-10053B		10 leads integrated; AHA, DB15M-Din3.0 			
EKG-10053C		10 leads integrated; IEC,DB15M-Banana4.0 			
EKG-10053D		10 leads integrated; AHA,DB15M-Banana4.0 			
EKG-10053E		10 leads integrated; IEC,DB15M-Clip end 			
EKG-10053F		10 leads integrated; AHA,DB15M-Clip end 			
EKG-10053G		10 leads integrated; IEC,DB15M-snap end 			
EKG-10053H		10 leads integrated; AHA,DB15M-snap end 			
Mortara Eli 250 cables					
Part No	OEM No	Description		Connector	Compatibility


EKG-10054A		10 leads integrated; IEC, DB15M-Din			
EKG-10054B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10054C		10 leads integrated; IEC,DB15M-Banana4.0			
EKG-10054D		10 leads integrated; AHA,DB15M-Banana4.0			
EKG-10054E		10 leads integrated; IEC,DB15M-Clip end			
EKG-10054F		10 leads integrated; AHA,DB15M-Clip end			
EKG-10054G		10 leads integrated; IEC,DB15M-snap end			
EKG-10054H		10 leads integrated; AHA,DB15M-snap end			
Dixtal EP-12 cables					
Part No	OEM No	Description		Connector	Compatibility
EKG-10055A		10 leads integrated; IEC, DB15M-Din			EP-12
EKG-10055B		10 leads integrated; AHA, DB15M-Din3.0			



RTARA
MENT: ELI
50/ELI250/
ELI210/ELI350

EKG-10055C		10 leads integrated; IEC,DB15M-Banana4.0			
EKG-10055D		10 leads integrated; AHA,DB15M-Banana4.0			
EKG-10055E		10 leads integrated; IEC,DB15M-Clip end			
EKG-10055F		10 leads integrated; AHA,DB15M-Clip end			
EKG-10055G		10 leads integrated; IEC,DB15M-snap end			
EKG-10055H		10 leads integrated; AHA,DB15M-snap end			
EKG-10056A		10 leads integrated; IEC, DB15M-Din			
EKG-10056B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10056C		10 leads integrated; IEC,DB15M-Banana4.0			
EKG-10056D		10 leads integrated; AHA,DB15M-Banana4.0			
EKG-10056E		10 leads integrated; IEC,DB15M-Clip end			Welch Allyn Cardioperfect Pro

EKG-10056F		10 leads integrated; AHA,DB15M-Clip end			
EKG-10056G		10 leads integrated; IEC,DB15M-snap end			
EKG-10056H		10 leads integrated; AHA,DB15M-snap end			

Welch Allyn CP100 CP200 EKG leads					
Part No	OEM No	Description		Connector	Compatibility
EKG-10057A		10 leads lead set; IEC, Din3.0			welch allyn CP100 CP200 ECG cable patient leads
EKG-10057B		10 leads lead set; AHA, Din3.0			
EKG-10057C		10 leads lead set; IEC, Banana4.0			
EKG-10057D	401129	10 leads lead set; AHA, Banana4.0			
EKG-10057E		10 leads lead set; IEC, clips			
EKG-10057F		10 leads lead set; AHA, clips			
EKG-10057G		10 leads lead set; IEC, snap			

EKG-10057H		10 leads lead set; AHA, snap			
EKG-10058A		10 leads lead set; IEC, Din3.0			
EKG-10058B		10 leads lead set; AHA, Din3.0			
EKG-10058C		10 leads lead set; IEC, Banana4.0			nova ecg leads
EKG-10058D		10 leads lead set; AHA, Banana4.0			
EKG-10058E		10 leads lead set; IEC, clips			
EKG-10058F		10 leads lead set; AHA, clips			
EKG-10058G		10 leads lead set; IEC, snap			
EKG-10058H		10 leads lead set; AHA, snap			
compatible Cardipia cables					
Part No	OEM No	Description		Connector	Compatibility
EKG-10059A		10 leads integrated; IEC, DB15M-Din			


EKG-10059B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10059C		10 leads integrated; IEC,DB15M-Banana4.0			 PanSW med Corp www.panswmed.com
EKG-10059D		10 leads integrated; AHA,DB15M-Banana4.0			
EKG-10059E		10 leads integrated; IEC,DB15M-Clip end			 PanSW med Corp www.panswmed.com
EKG-10059F		10 leads integrated; AHA,DB15M-Clip end			
EKG-10059G		10 leads integrated; IEC,DB15M-snap end			
EKG-10059H		10 leads integrated; AHA,DB15M-snap end			

compatible Bruker cables					
Part No	OEM No	Description		Connector	Compatibility
EKG-10060A		10 leads integrated; IEC, DB15M-Din			
EKG-10060B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10060C		10 leads integrated; IEC,DB15M-Banana4.0			

EKG-10060D		10 leads integrated; AHA,DB15M-Banana4.0		
EKG-10060E		10 leads integrated; IEC,DB15M-Clip end		
EKG-10060F		10 leads integrated; AHA,DB15M-Clip end		
EKG-10060G		10 leads integrated; IEC,DB15M-snap end		
EKG-10060H		10 leads integrated; AHA,DB15M-snap end		


compatible Zoncare cables					
Part No	OEM No	Description		Connector	Compatibility
EKG-10061A		10 leads integrated; IEC, DB15M-Din			
EKG-10061B		10 leads integrated; AHA, DB15M-Din3.0			
EKG-10061C		10 leads integrated; IEC,DB15M-Banana4.0			
EKG-10061D		10 leads integrated; AHA,DB15M-Banana4.0			
EKG-10061E		10 leads integrated; IEC,DB15M-Clip end			


EKG-10061F		10 leads integrated; AHA,DB15M-Clip end			
EKG-10061G		10 leads integrated; IEC,DB15M-snap end			
EKG-10061H		10 leads integrated; AHA,DB15M-snap end			
compatible Conmed ECG leads					
Part No	OEM No	Description		Connector	Compatibility
EKG-10062A		10 leads; IEC, DB15M-Din3.0			
EKG-10062B		10 leads; AHA, DB15M-Din3.0			
EKG-10062C		10 leads; IEC,DB15M-Banana4.0			
EKG-10062D		10 leads; AHA,DB15M-Banana4.0			
EKG-10062E		10 leads; IEC,DB15M-Clip end			
EKG-10062F		10 leads; AHA,DB15M-Clip end			
EKG-10062G		10 leads; IEC,DB15M-snap end			

EKG-10062H		10 leads; AHA,DB15M-snap end			
compatible CU ECG cable					
Part No	OEM No	Description		Connector	Compatibility
EKG-10063A		10 leads; IEC, DB15M-Din3.0			CU medical
EKG-10063B		10 leads; AHA, DB15M-Din3.0			
EKG-10063C		10 leads; IEC,DB15M-Banana4.0			
EKG-10063D		10 leads; AHA,DB15M-Banana4.0			
EKG-10063E		10 leads; IEC,DB15M-Clip end			
EKG-10063F		10 leads; AHA,DB15M-Clip end			
EKG-10063G		10 leads; IEC,DB15M-snap end			
EKG-10063H		10 leads; AHA,DB15M-snap end			
compatible Sadaat ECG cable/ Heyer ECG cable					
Part No	OEM No	Description		Connector	Compatibility

EKG-10064A		10 leads; IEC, DB15M-Din3.0			compatible with Sadaat medical
EKG-10064B		10 leads; AHA, DB15M-Din3.0			
EKG-10064C		10 leads; IEC, DB15M-Banana4.0			
EKG-10064D		10 leads; AHA, DB15M-Banana4.0			
EKG-10064E		10 leads; IEC, DB15M-Clip end			
EKG-10064F		10 leads; AHA, DB15M-Clip end			
EKG-10064G		10 leads; IEC, DB15M-snap end			
EKG-10064H		10 leads; AHA, DB15M-snap end			
compatible GE CAM 14 main cable for acquisition box					
CAM-TRK-14		original part number: 2016560-001 (short); 2016560-003 (long wires)			

ADAPTERS

EKG-ADP-1	Adaptor banana 4mm - Snap			Adaptor banana 4mm - Snap with wires	
EKG-ADP-2	Adaptor banana 4mm- Clip			Adaptor banana 4mm-Clip	
EKG-ADP-3	Adaptor banana 4mm- SNAP			Adaptor banana 4mm-SNAP	
EKG-ADP-4	4.0 banana to crocodile clip		4.0 banana to crocodile clip		
EKG-ADP-5	 	 	Universal Flat Clip Adapter w Bannana ECG Electrodes; Philips part number 989803166031	4.0 banana to button electrodes, tap electrodes	
EKG-ADP-6	 	 		4.0 banana to button electrodes, tap electrodes	
EKG-ADP-7	 	 	NKD-H/ADAPTER2 GE 2104729-001	4.0 banana to button electrodes, tap electrodes	
EKG-ADP-8	Adaptor banana 4mm - Snap			Adaptor banana 4mm - Snap with wires	

EKG-ADP-9	Adaptor 2.0 banana to 4mm-SNAP	 	Adaptor 2.0 banana 4mm-	
EKG-ADP-10	3.0 DIN to crocodile clip	 	3.0 DIN to crocodile clip	