


Email: sales@panswmed.com
 Tel: 86-20-28961025
 Skype: [panmedical](https://www.skype.com/partners/panmedical)


PANSW NIBP Connectors


Part No	Description	Picture	Compatibility	Original part number
CNT-01	plastic		all kinds of disposable nibp cuffs; Male Slip Luer to 5/32 in. ID Tube, Plastic	300666
CNT-02	plastic		Nihon Kohden dual tube system interconnect hoses; other dual tube systems	GE Part # 300829
CNT-03 (F) CNT-03 (M)	    Datex Ohmeda submin connectors Mated Pair Submin to 1/8" ID Barbed, Plastic Compatibility: Datex-Ohmeda Part # 330058 Artema NIBP tubing system Female Submin to 1 /8" ID Barbed, Plastic Compatibility: GE, Marquette Protocol Part # 330064 Male Submin to 1 /8" ID Barbed, Plastic Compatibility: GE, Marquette, Datex-Ohmeda Part # 330091			
CNT-04			Fukuda Kontron, Schiller	Use with CNT-06

CNT-05	neonate NIBP connector for all the marks neonate system		5/32" ID Barbed to Female Slip Luer , Metal Compatibility: Part # 300669 Use with CNT-18	
CNT-06				Use with CNT-04
CNT-07(F) CNT-07(M)			Compatibility: Philips/HP,Sieme ns, Datascope, Colin Part # 330060 Compatibility: Philips/HP,Sieme ns, Datascope, Colin; Mindray Part # 330059 Spacelabs 714- 0045-00;	5/32" ID Barbed to Female Bayonet, Metal Male Bayonet to 5/32" ID Barbed, Metal Spacelabs: 714-0038- 00
CNT07FF			outer diameter: 5.3mm	
CNT-08				Nihon Kohden adult & pediatric

CNT-08N			Nihon Kohden neonate
CNT-09			GE Marquette, Ohmeda, Welch Allyn Part# 300619 1/8" ID Barbed to Female Screw, Metal; Compatibility: DINAMAP®; Use with CNT-19/CNT-11
CNT-10			neonate NIBP connector for disposable NIBP cuff
CNT-11			Compatibility: DINAMAP® Part # 300664
CNT-12			Compatibility: DINAMAP® Pro, Compact, MPS Series Monitors
			Male Screw to 5/ 32" ID Barbed, Plastic ; Use with CNT-09
			330073

CNT-13			Compatibility: Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120	
CNT-14			GE P/N:300670	GE P/N: 300670 size: 4mm(5/32 in.) to 2.5mm Two end thread
CNT-15				use with CNT-20
CNT-16				
CNT-17				

CNT-18		5/32" ID Barbed to Male Slip Luer, Metal Part # 300667 Use with CNT-05
CNT-19		Metal version of CNT-11; Use with CNT-09
CNT-20		Use with CNT-15
CNT-21		Compatibility: Spacelabs, Datascope, Colin; Female Locking Luer to 5/32" ID Barbed, Plastic; Part# 300668; Use with CNT-22
CNT-22		Compatibility: Spacelabs, Datascope, Colin; Part # 300665 5/32" ID Barbed to Male Locking Luer, Plastic; Use with CNT-21
CNT-23		GE; Datex Ohmeda NIBP connector (monitor side)

CNT-24			plastic version of CNT-07(M); Use with CNT-07(F) barb OD: 5mm
CNT-24N			for neonate use, barb OD: 4.1mm
CNT-25			5/32" ID Barbed to Male Slip Luer, plastic Use with CNT-05/CNT-26
CNT-26			Female slip luer connector; plastic; Use with CNT- 25/CNT-18
CNT-27			

CNT-28			Pace Tech; marquette
CNT-29	 		Creative, Neusoft, Jeteem; Kernel
CNT-30A			nibp adapter dual tubes to single tube
CNT-30B			nibp adapter dual tubes to single tube
CNT-30C			nibp adapter dual tubes to single tube


CNT-31			dual to single nibp hose adatper
CNT-32			connector for monitor side
CNT-33			
CNT-34			omron blood meter
CNT-35			philips machine connector

CNT-36			Drager machine side connector
CNT-37A			CPC style connector, for machine port replacement use
CNT-37B			CPC style connector, for hose replacement use
CNT-38			
CNT-39			drager nibp cuff connector

CNT-51			drager nibp hose connector
CNT-40			hose connector for philips new
CNT-41			hose connector for philips new
CNT-42			welch allyn nibp hose connector for machine side
CNT-43			
CNT-44			new GE neonate dual hoses nibp connector

CNT-45			welch allyn new connector dual hoses
CNT-46			welch allyn new connector single hoses
CNT-47			Mindray NIBP connector, metal one to use on new beneview and other models.
CNT-48	 PanSW Med Corp www.panswmed.com		nihon kohden opv nibp hose connector, adult/child type
CNT-49			nihon kohden opv nibp hose connector, neonate type
CNT-50	 PanSW Med Corp www.panswmed.com		machine connector for infinium, newtech...
CNT-51	 PanSW Med Corp www.panswmed.com		drager nibp hose connector

CNT-52		dual hose nibp connector
CNT-53		omron nibp connector bent, connector to machine OD:5.0mm
CNT-54		omron nibp connector bent, connector to machine OD:6.0mm
CNT-55		male connector
CNT-56		female connector

CNT-57			single hose use, to work with CNT-43
CNT-58			DUAL hose use, to work with CNT-43